

# Introducing the NEW ACCURATE ONE STEP URINE CUP



Now Testing  
for Synthetic  
Marijuana!

## **EXCLUSIVE** to LaneWorkSafe

The answer to Australia's workplace drug testing needs

### **Because keeping staff safe is the priority...**

The Accurate One Step Urine Cup was designed out of the increasing concern and frustration safety professionals faced in finding a solution for the on-site detection of Synthetic Marijuana, commonly known as "Kronic".

The all-in-one innovative cup design, unique to the Accurate One Step Urine Cup, is the necessary solution for workplace safety programs allowing for the reliable on-site detection of 7 major drug groups.

The device is specially designed to commence testing once the donor voids into the cup. This means less time is spent on waiting for results, so workers can get back quickly and safely to the job at hand.



This device has been issued with a Certificate of Compliance by an Australian NATA Accredited Laboratory proving it meets AS/NZ 4308:2008 cut-off levels as required.

## THE FACTS: DID YOU KNOW?

- Synthetic Marijuana is one of the latest psychoactive substances distributed in Australia
- JWH-018 and JWH-073 are the primary synthetic cannabinoids responsible for the euphoric and psychoactive effects that imitate marijuana
- Until recently, Synthetic Marijuana was unable to be detected in traditional immunoassay urine screening devices
- According to the latest Australian Crime Commission figures, seizures of synthetic drugs have increased by 64.3% from 2013 to 2014\*  
\*[http://www.druginfo.adf.org.au/attachments/article/1363/FS\\_Workplace\\_Sep2015.pdf](http://www.druginfo.adf.org.au/attachments/article/1363/FS_Workplace_Sep2015.pdf)

### Use with Confidence

- ✓ Quality assured
- ✓ AS/NZ 4308:2008 Compliant
- ✓ Results in 2 - 3 minutes
- ✓ Cost Effective
- ✓ Easy to use
- ✓ Easy to interpret
- ✓ No leakage

## DOCUMENTED RELIABILITY AND ACCURACY

There is well documented evidence that urine drug screening remains as one of the most reliable methods of on-site drug screening in Australia. According to the National Association of Testing Authorities (NATA) there are no clearly defined cut-off concentrations for devices published in AS4760:2007 as there are for urine devices in AS/NZS 4308:2008.\* The New Accurate One Step Urine Cup has the capabilities of testing for additional tailored substances which can include but are not limited to; Oxycodone, Buprenorphine, Methadone, Fentanyl, Tramadol and Ketamine.

\*<http://laneworksafe.com.au/nata-report-important-findings-saliva-accreditation/>

Types of Drugs	Abbreviations	Cut Off Levels
Synthetic Marijuana	K2	30 ng/ml
Amphetamine	AMP	300 ng/mL
Methamphetamine	mAMP	300 ng/mL
Cocaine	COC	300 ng/mL
Marijuana	THC	50 ng/mL
Benzodiazepines	BZO	200 ng/mL
Morphine	MOP	300 ng/mL

The Accurate One Step Urine Cup is compliant with AS/NZ 4308:2008 Cut-Off Levels, as per the diagram.

## Also Available Confirmatory Testing Kits

Got a non-negative result? Accurate One Step Urine Cup can provide evacuator kits for laboratory GCMS testing. Confirmatory testing Kits can be purchased separately. The kit includes 3 pieces: 1 cup, 1 lid and vials for collection.

**TO GET IN TOUCH CALL**  
**1800 429 219** or email  
**[sales@laneworksafe.com.au](mailto:sales@laneworksafe.com.au)**

**YOU CAN ALSO FIND US ON**  
[laneworksafe.com.au](http://laneworksafe.com.au)  
[linkedin.com/company/laneworksafe-pty-ltd](https://www.linkedin.com/company/laneworksafe-pty-ltd)  
[twitter.com/LaneWorkSafe](https://twitter.com/LaneWorkSafe)

