

# Introducing the **ACCURATE ORAL FLUID DRUG SCREEN 9-PANEL**



NOW TESTING FOR  
**SYNTHETIC  
MARIJUANA  
& ALCOHOL**

## **EXCLUSIVE** to LaneWorkSafe

The answer to Australia's workplace drug testing needs.

### *Accurate, Reliable, Cost Effective.*

The Accurate Oral Fluid Drug Screen 9-Panel is one of the most comprehensive saliva-based drug tests available on the market, screening for a total of nine drug groups onsite. Its range of detection includes Synthetic Marijuana, Oxycodone, and Alcohol, reinforcing safety programs to uphold the highest standards and latest technologies in onsite testing.

LaneWorkSafe addresses the option for workplace use of saliva or oral fluid drug testing with an integrated, user friendly, and compliant drug test. Equipped with an enhanced mint-flavored sponge to facilitate plentiful sample collection, this device will allow for expanded testing within workplaces across Australia.



This device has been issued with a Compliance Certificate proving it meets AS/NZ 4760:2019 cut off levels and requirements.

## THE FACTS: DID YOU KNOW?

- Recent research has estimated that 2.5 million days are lost annually due to drinking and drug use.
- Alcohol and other drugs cost Australian workplaces an estimated \$6 billion per year in lost productivity.
- 1 in 10 workers say they have been affected by a co-worker's use of alcohol.
- Use of Synthetic Cannabinoids has been associated with a range of adverse health effects. Including; nausea, paranoia, seizures, and in some cases death.

### Use with Confidence

- ✓ Quality assured
- ✓ AS/NZ 4760:2019 Compliant
- ✓ Easy to use
- ✓ Easy to interpret
- ✓ Cost-effective
- ✓ Tamper proof
- ✓ Results in minutes

## DOCUMENTED RELIABILITY AND ACCURACY

The LaneWorkSafe Accurate Oral Fluid Drug Screen has been certified by a NATA Accredited Laboratory as meeting the requirements of AS/NZ 4760:2019. This certification ensures the saliva device meets with Australian Standards demonstrating its reliability and suitability for use in workplaces.

Types of Drugs	Abbreviations	Cut-Off Levels
Synthetic Marijuana	K2	25 ng/mL
Amphetamine	AMP	50 ng/mL
Methamphetamine	mAMP	50 ng/mL
Cocaine	COC	50 ng/mL
Marijuana	THC	15 ng/mL
Benzodiazepines	BZO	50 ng/mL
Morphine	MOP	50 ng/mL
Alcohol	ALC	0.2% BAC
Oxycodone	OXY	40 ng/mL

### TO GET IN TOUCH

Call **1800 429 219**  
or email  
**[sales@laneworksafe.com.au](mailto:sales@laneworksafe.com.au)**