



LaneWorkSafe

Worksafe – Workwell

Split-Specimen Cup™

Innovative drug screen cup design eliminates urine handling and donor tampering, providing a unique and efficient approach for on-site drug screening.



*This product is compliant with AS4308:2008 cut-off levels.
Compliance Certificate available on request.*

Use with confidence

- Split-specimen collection
- No urine handling - ready for GC/MS
- Screw cap with audible lock
- Female donor friendly
- Screens for 6 drugs and 6 adulterants
- Built-in temperature strip
- Easy to photocopy results
- Results in around 2 minutes

**Used by NSW Police
for officer testing.**



LaneWorkSafe
Worksafe – Workwell

LaneWorkSafe Pty Ltd | ACN 113 839 713 | PO Box 266 Westgate NSW 2048
T 1800 429 219 F 02 9559 1559 www.laneworksafe.com.au

Point of care medical diagnostic screening products